



OVERVIEW

The VUW-7 Digital Ultrasound Wattmeter provides quick and accurate verifying and calibrating readout of therapeutic ultrasound devices in the field or hospital settings.

WHAT MAKES THE UW-7 UNIQUE?

- Heavy cast frame firmly supports rugged straight bar load cell. Provides for faster technician on-site visits and protects the load from cell damage during transit and handling.
- Support Ring Sets provide precise sound head centering and alignment for accurate

★ **Compact** ★ **Easy to Use** ★ **Best Value**

PRODUCT HIGHLIGHTS

- Special temperature compensated load cell design
- High Accuracy 0.05%
- 0.10 watt to 30 watt power measurement range
- Fast auto zeroing button
- Extremely fast 1.5 second readout stabilization
- Measures pulsed or continuous power
- 2.65" (6.7 cm) diameter sound head maximum
- High contrast back lit display
- Selectable readout in watts or grams
- 120 or 240 VAC 50-60 Hz standard
- 100 mW resolution
- Accuracy = $\pm 3\%$ + one count reading

SPECIFICATIONS

Unit Dimensions 8.5" W, 12" D, 7.25" H

Unit weight 8.25 LBS (3.74 kg)

Selectable readouts Watts or Grams

Function Check Weights 75g or 200g

Ordering Information

| Part No: | Description |
|----------|--------------------------|
| VUW-7 | New Ultrasound Wattmeter |

Standard Accessories:

| | |
|----------|--------------------|
| VU7-CASE | Soft Carrying Case |
|----------|--------------------|

Optional Accessories:

- Optional RS-232 cable available
- Optional USB cable available
- Optional high capacity lithium-ion battery power module available



For other popular Netech products, visit www.NetechCorp.US