

## Netech Defibrillator Analyzer Models Comparison

	<b>DELTA 1500</b>	<b>DELTA 3000</b>	<b>DELTA 3300 (New)</b>
<b>TEST FUNCTIONS</b>	<ul style="list-style-type: none"> <li>Tests all AEDs (Automated External Defibrillators) including pulsed AEDs.</li> </ul>	<ul style="list-style-type: none"> <li>Tests all AEDs</li> <li>Tests all monophasic and biphasic defibrillators and transcutaneous Pacemakers.</li> <li>Tests all External Pacemakers.</li> </ul>	<ul style="list-style-type: none"> <li>Tests all AEDs</li> <li>Tests all monophasic, biphasic and pulsed defibrillators and transcutaneous Pacemakers.</li> </ul>
<b>MAIN FEATURES</b>	<ul style="list-style-type: none"> <li>Compact</li> <li>Battery operated</li> <li>Easy to use</li> </ul>	<ul style="list-style-type: none"> <li>Compact, lightweight, and rugged</li> <li>Simultaneous display of all test measurements</li> <li>Simple menu driven operation Energy measurement</li> <li>Pacemaker analyzer (Transcutaneous/ Atrial and Ventricular).</li> <li>Demand sensitivity, Immunity test, refractory period, waveform width</li> <li>ECG arrhythmia simulator</li> </ul>	<ul style="list-style-type: none"> <li>Compact and lightweight (3 lbs)</li> <li>Simple and easy to use</li> <li>Li-ion rechargeable battery</li> <li>Large LCD display with backlight</li> <li>Incredible value</li> <li>Energy measurement; sync time, charge time, waveform storage/playback</li> <li>Transcutaneous pacemaker analyzer; Demand sensitivity, Immunity test, refractory period, waveform width.</li> </ul>
<b>ECG SIMULATOR</b>	<ul style="list-style-type: none"> <li>Simulates normal sinus rhythm, ventricular fibrillation, and ventricular tachycardia waveforms</li> </ul>	<ul style="list-style-type: none"> <li>12 lead normal sinus rhythm and pulse</li> <li>Clinical and performance waveforms</li> </ul>	<ul style="list-style-type: none"> <li>ECG arrhythmia simulator</li> <li>12 lead normal sinus rhythm, performance waveform and pacer pulse.</li> <li>Clinical and performance waveforms.</li> </ul>
<b>SERIAL PORT/ TEST DATA SAVE CAPABILITY</b>	None	<ul style="list-style-type: none"> <li>Serial port</li> <li>Dedicated print key for immediate output of test data.</li> </ul>	<ul style="list-style-type: none"> <li>Serial port</li> <li>Capable of saving up to 5,000 asset IDs and Test records.</li> <li>Upload to PC for printing/record keeping.</li> <li>Free ProCom software included</li> </ul>
<b>DIMENSIONS</b>	7.8" X4.3" X 2" (19.7 X10.8 X5 Cm)	10" x 11" x 2.25" (25.4 x 27.9 x 5.7 cm)	10.25 x 6.25 x 2.5 Inch. (26 X 16 x 6.4 cm)
<b>WEIGHT</b>		4.0 lbs (1.8 kg)	3.3 lbs (1.5 kg)
<b>PRICE</b>	<del>\$1139.00</del>	<del>\$3995</del> (Clearance)	<del>\$4150</del> <b>Special Introductory</b>
<i>Promotional Price</i>	<i>Call</i> 1800-547-6557 Email <a href="mailto:Sales@Netech.org">Sales@Netech.org</a>	<i>Call</i> 1800-547-6557 Email <a href="mailto:Sales@Netech.org">Sales@Netech.org</a>	<i>Call</i> 1800-547-6557 Email <a href="mailto:Sales@Netech.org">Sales@Netech.org</a>